

Santé publique Ontario

Public Health Laboratory

Laboratoire de santé publique

Labstract - March 2014

Acid Fast Bacilli (AFB) Smear Reports - Correlation with the Canadian Tuberculosis Standards, 7th Edition

To Health Care Providers:

The Public Health Ontario Laboratories (PHOL) AFB smear reporting protocol reflects an internationally accepted semi-quantitative terminology and protocol for enumerating AFB observed in smears. The following chart compares AFB smear results from PHOL enumeration methods to those published in the Canadian Tuberculosis Standards (CTS) -7^{th} Edition:

Number of AFB seen in smear	PHOL Reporting ²	CTS Reporting
None	No AFB seen	No AFB seen
3-9 / smear	Few ¹	Inconclusive, repeat ¹
1-9 / ten fields	1+	1+ (rare) ³
1-9 / field	2+	2+ (few)
10-90 / field	3+	3+ (moderate)
> 90 / field	4+	4+ (numerous)

¹Few as designated by PHOL is reported only after a repeat second smear is read confirming the presence of acid-fast bacilli on the smear. PHOL does not use the "inconclusive, repeat" terminology.

Reference: Canadian Tuberculosis Standards – 7th Edition, Canadian Respiratory Journal Volume 20, Supplement A July/August 2013



²The criteria set in the Canadian Tuberculosis Standards are enumerated by reading smears at 250x magnification. PHOL smears are read at 200x magnification with the same number of fields as in the CTS. Results correlate with the same enumeration as indicated in the chart above.

³The PHOL does not use the terminology as described in the CTS (rare, few etc). Please refer to the chart above regarding smear enumeration.

Acid Fast Bacilli Reports (Continued)

For further information:

- Contact the PHOL Customer Service Centre at 416-235-6556 or 1-877-604-4567 (toll-free), or by email at customerServiceCentre@oahpp.ca
- For PHOL specimen collection information and previous Labstracts, refer to http://www.publichealthontario.ca/Labs
- The current version of the PHOL General Test Requisition and other forms are available at http://www.publichealthontario.ca/Requisitions
- To subscribe to future Labstracts, email labstracts@oahpp.ca
- To register for Autofax and receive laboratory reports by fax directly from our laboratory information system as soon as they are released, contact the PHOL Customer Service Centre.